



TrioDocs

Version: 0.4.0

Date: August 23, 2025

Download the latest version on:  
<https://triodocs.org>

# Overrides

# Overrides

## What is an Override?

An Override is a set of parameters you can place on Trio temporarily. With an Override, you can alter a variety of variables to increase or decrease insulin delivery:

- Basal Rate Adjustment
- Inversely change ISF, CR, or both
- Override Target
- Disable SMBs
- Override Max SMB Minutes
- Set a time frame for the override or enable indefinitely until you cancel

There are multiple cases when an Override may be useful. Below you will find a few examples. This list is not all-inclusive and you may find additional uses for Overrides as you use Trio.

In addition to Overrides, Trio also has a more automated, less customizable feature called [Temporary Targets](#).

---

## How Do I Set an Override?

Overrides can be found under the Adjustments tab at the bottom of the Trio app screen.

You can also find your saved Overrides listed [here](#).

Overrides can be set remotely using [Loop Follow](#) or automated using [Shortcuts](#).

---

## Examples of Use