



TrioDocs

Version: 0.2.1

Date: May 14, 2025

Download the latest version on:  
<https://docs.diy-trio.org>

# Notifications

# Notifications

Allows configuration of Trio notifications. Adjust based on your preferences.

Notifications can be enabled/disabled by toggling `Allow Notifications` on/off in `iPhone Settings > Notifications > Trio`.

## Glucose

### Show glucose on the app badge

This will add your current glucose on top of your Trio icon.

### Always Notify Glucose

A notification will be triggered every time your glucose is updated in Trio.

### Always play alert sound

This will cause a sound to be triggered by every Trio notification.

### Also add source info

The source of the glucose reading will be added to the notification.

## Low

Set this to the highest value you want a low glucose notification to be triggered at.

## High

Set this to the lowest value you want a high glucose notification to be triggered at.

## Carbs Required Threshold

Trio may predict a low blood sugar event and recommend consuming a suggested amount of carbs. `Carbs Required Threshold` determines the minimum amount of carbs required before Trio sends a notification.