



TrioDocs

Version: 0.2.1

Date: May 15, 2025

Download the latest version on:  
<https://triodocs.org>

# Migration from iAPS

# Coming from iAPS

When you move your settings from iAPS to Trio, some of them might look different even if they have the same or similar names. That's because we changed some of the numbers from decimals to percentages to make them easier to understand.

 **Example**

In iAPS you might set "Autosens Maximum" to 1.2. In Trio, we show "Autosens Max" as 120%, which means the same thing, but is easier to understand.

Some setting names have also changed. We did this to make things less confusing. To help you out, we added definitions right in the app. Just tap the question mark icon next to any setting to see what it means. Healthcare Professionals and users can also find those settings and explanations [in this section of the docs](#).

In iAPS, some settings had limits (called guardrails) that weren't always visible, and others had no limits at all. In Trio, we've made those guardrails easier to see and added a few more where needed, to help prevent settings that could lead to unsafe outcomes.

The charts below will help you see which setting names or formats have changed, and where to find them in both iAPS and Trio.

 **Tip**

To convert a value from a decimal to a percentage, multiply it by 100.  
To convert from a percentage to a decimal, divide by 100.

## Prepare Trio

| Trio Name     | Setting Format (example)    | Location in App                       | iAPS Name      | Setting Format (example)    | Location in App                    |
|---------------|-----------------------------|---------------------------------------|----------------|-----------------------------|------------------------------------|
| Glucose Units | selection (mg/dL or mmol/L) | Settings → Therapy → Units and Limits | Target Glucose | selection (mg/dL or mmol/L) | Settings → OpenAPS → Glucose units |

## Therapy Settings

| Trio Name | Setting Format (example) | Location in App | iAPS Name | Setting Format (example) | Location in App |
|-----------|--------------------------|-----------------|-----------|--------------------------|-----------------|
|-----------|--------------------------|-----------------|-----------|--------------------------|-----------------|

|                              |                                        |                                            |                              |                                        |                                  |
|------------------------------|----------------------------------------|--------------------------------------------|------------------------------|----------------------------------------|----------------------------------|
| <b>Glucose Targets</b>       | decimal<br>(100 mg/dL   5.5 mmol/L)    | Settings → Therapy → Glucose Targets       | <b>Target Glucose</b>        | decimal<br>(100 mg/dL   5.5 mmol/L)    | Settings → Target Glucose        |
| <b>Basal Rates</b>           | decimal<br>(1.0 U/hr)                  | Settings → Therapy → Basal Rates           | <b>Basal Profile</b>         | decimal<br>(1.0 U/hr)                  | Settings → Basal Profile         |
| <b>Carb Ratios</b>           | decimal<br>(10 g/U)                    | Settings → Therapy → Carb Ratios           | <b>Carb Ratios</b>           | decimal<br>(10 g/U)                    | Settings → Carb Ratios           |
| <b>Insulin Sensitivities</b> | decimal<br>(54 mg/dL/U   3.0 mmol/L/U) | Settings → Therapy → Insulin Sensitivities | <b>Insulin Sensitivities</b> | decimal<br>(54 mg/dL/U   3.0 mmol/L/U) | Settings → Insulin Sensitivities |

## Delivery Limits

| Trio Name                       | Setting Format (example) | Location in App                                          | iAPS Name                | Setting Format (example) | Location in App                                     |
|---------------------------------|--------------------------|----------------------------------------------------------|--------------------------|--------------------------|-----------------------------------------------------|
| <b>Max IOB</b><br>.....         | decimal<br>(2 U)         | Settings → Therapy → Units and Limits → Max IOB<br>..... | <b>Max IOB</b><br>.....  | decimal<br>(2)           | Settings → OpenAPS → Max IOB<br>.....               |
| <b>Max Bolus</b>                | decimal<br>(10 U)        | Settings → Therapy → Units and Limits → Max Bolus        | <b>Max Bolus</b>         | decimal<br>(10)          | Settings → Pump Settings → Max Bolus                |
| <b>**Max Basal Rate</b>         | decimal<br>(2 U/hr)      | Settings → Therapy → Units and Limits → Max Basal Rate   | <b>Max Basal</b>         | decimal<br>(2)           | Settings → Pump Settings → Max Bolus                |
| <b>Max COB</b><br>.....         | decimal<br>(120 g)       | Settings → Therapy → Units and Limits → Max COB<br>..... | <b>Max COB</b><br>.....  | decimal<br>(120)         | Settings → OpenAPS → Max COB<br>.....               |
| <b>Minimum Safety Threshold</b> | decimal<br>(60 mg/dL)    | Settings → Therapy → Units and Limits → Max IOB<br>..... | <b>Threshold Setting</b> | decimal<br>(60)          | Settings → Dynamic ISF → Threshold Setting<br>..... |

# Algorithm Settings

## Autosens

| Trio Name           | Setting Format<br>(example) | Location in App                                         | iAPS Name               | Setting Format<br>(example) | Location in App                                |
|---------------------|-----------------------------|---------------------------------------------------------|-------------------------|-----------------------------|------------------------------------------------|
| <b>Autosens Min</b> | percentage<br>(70%)         | Settings →<br>Algorithm →<br>Autosens →<br>Autosens Min | <b>Autosens Minimum</b> | decimal<br>(0.7)            | Settings →<br>OpenAPS →<br>Autosens<br>Minimum |
| <b>Autosens Max</b> | percentage<br>(120%)        | Settings →<br>Algorithm →<br>Autosens →<br>Autosens Max | <b>Autosens Maximum</b> | decimal<br>(1.2)            | Settings →<br>OpenAPS →<br>Autosens<br>Maximum |

## SMB (Super Micro Bolus)

| Trio Name                              | Setting Format<br>(example) | Location in App                                                                             | iAPS Name                             | Setting Format<br>(example) | Location in App                                                 |
|----------------------------------------|-----------------------------|---------------------------------------------------------------------------------------------|---------------------------------------|-----------------------------|-----------------------------------------------------------------|
| <b>Enable SMB Always</b>               | toggle<br>(On/Off)          | Settings →<br>Algorithm → Super<br>Micro Bolus (SMB)<br>→ Enable SMB<br>Always              | <b>Enable SMB Always</b>              | toggle<br>(On/Off)          | Settings →<br>OpenAPS →<br>Enable SMB<br>Always                 |
| <b>Allow SMB with High Temp Target</b> | toggle<br>(On/Off)          | Settings →<br>Algorithm → Super<br>Micro Bolus (SMB)<br>→ Allow SMB With<br>High Temptarget | <b>Allow SMB With High Temptarget</b> | toggle<br>(On/Off)          | Settings →<br>OpenAPS →<br>Allow SMB With<br>High<br>Temptarget |
| <b>Enable UAM</b>                      | toggle<br>(On/Off)          | Settings →<br>Algorithm → Super<br>Micro Bolus (SMB)<br>→ Enable UAM                        | <b>Enable UAM</b>                     | toggle<br>(On/Off)          | Settings →<br>OpenAPS →<br>Enable UAM                           |

|                                                 |                     |                                                                                            |                                           |                  |                                                          |
|-------------------------------------------------|---------------------|--------------------------------------------------------------------------------------------|-------------------------------------------|------------------|----------------------------------------------------------|
|                                                 |                     | Minutes                                                                                    | Enable UAM                                | decimal<br>(30)  | Settings →<br>OpenAPS → Max<br>SMB Basal<br>Minutes      |
| <b>Max UAM<br/>Basal<br/>Minutes</b>            | decimal<br>(30 min) | Settings → Algorithm →<br>Super Micro Bolus<br>(SMB) → Max UAM<br>Basal Minutes            | <b>Max UAM<br/>SMB Basal<br/>Minutes</b>  | decimal<br>(30)  | Settings →<br>OpenAPS → Max<br>UAM SMB Basal<br>Minutes  |
| <b>Max Allowed<br/>Glucose Rise<br/>for SMB</b> | percentage<br>(20%) | Settings → Algorithm →<br>Super Micro Bolus<br>(SMB) → Max Allowed<br>Glucose Rise for SMB | <b>Max Delta-BG<br/>Threshold<br/>SMB</b> | decimal<br>(0.2) | Settings →<br>OpenAPS → Max<br>Delta-BG<br>Threshold SMB |

## Target Behavior

| Trio Name                                          | Setting<br>Format<br>(example) | Location in App                                                                          | iAPS Name                                             | Setting<br>Format<br>(example) | Location in App                                                     |
|----------------------------------------------------|--------------------------------|------------------------------------------------------------------------------------------|-------------------------------------------------------|--------------------------------|---------------------------------------------------------------------|
| <b>High Temp<br/>Target Raises<br/>Sensitivity</b> | toggle<br>(On/Off)             | Settings →<br>Algorithm →<br>Target Behavior →<br>High Temp Target<br>Raises Sensitivity | <b>High<br/>Temptarget<br/>Raises<br/>Sensitivity</b> | toggle<br>(On/Off)             | Settings →<br>OpenAPS →<br>High Temptarget<br>Raises<br>Sensitivity |
| <b>Low Temp<br/>Target Lowers<br/>Sensitivity</b>  | toggle<br>(On/Off)             | Settings →<br>Algorithm →<br>Target Behavior →<br>Low Temp Target<br>Lowers Sensitivity  | <b>Low<br/>Temptarget<br/>Lowers<br/>Sensitivity</b>  | toggle<br>(On/Off)             | Settings →<br>OpenAPS →<br>Low Temptarget<br>Lowers<br>Sensitivity  |
| <b>Sensitivity<br/>Raises Target</b>               | toggle<br>(On/Off)             | Settings →<br>Algorithm →<br>Target Behavior →<br>Sensitivity Raises<br>Target           | <b>Sensitivity<br/>Raises Target</b>                  | toggle<br>(On/Off)             | Settings →<br>OpenAPS →<br>Sensitivity<br>Raises Target             |
| <b>Resistance<br/>Lowers Target</b>                | toggle<br>(On/Off)             | Settings →<br>Algorithm →<br>Target Behavior →<br>Resistance Lowers<br>Target            | <b>Resistance<br/>Lowers Target</b>                   | toggle<br>(On/Off)             | Settings →<br>OpenAPS →<br>Resistance<br>Lowers Target              |
|                                                    |                                |                                                                                          |                                                       |                                |                                                                     |

|                                  |                                           |                                                                              |                                  |                                           |                                                          |
|----------------------------------|-------------------------------------------|------------------------------------------------------------------------------|----------------------------------|-------------------------------------------|----------------------------------------------------------|
| Half Basal<br>Exercise<br>Target | decimal<br>(160 mg/dL<br>  8.9<br>mmol/L) | Settings → Algorithm<br>→ Target Behavior →<br>Half Basal Exercise<br>Target | Half Basal<br>Exercise<br>Target | decimal<br>(160 mg/dL<br>  8.9<br>mmol/L) | Settings →<br>OpenAPS → Half<br>Basal Exercise<br>Target |
|----------------------------------|-------------------------------------------|------------------------------------------------------------------------------|----------------------------------|-------------------------------------------|----------------------------------------------------------|