



TrioDocs

Version: 0.4.0

Date: March 24, 2026

Download the latest version on:
<https://triodocs.org>

Compatible Devices

Compatible Devices

Below is the list of compatible devices available for use with Trio.

Phone

Trio is compatible with phones running on **iOS 17 or newer**.

These iPhones are compatible with **iOS version 17 and newer**:

- iPhone 17, all variants
- iPhone 16*, all variants
- iPhone 15, all variants
- iPhone 14, all variants
- iPhone 13, all variants
- iPhone 12, all variants
- iPhone 11, all variants
- iPhone XS Max
- iPhone XS
- iPhone XR
- iPhone SE (3rd generation; 2022 first release)
- iPhone SE (2nd generation; 2020 first release)

iPhone 16 with Atlas/InPlay DASH pods

Using an iPhone 16 with Atlas/InPlay DASH pods can cause communication issues, but SAW/TWI DASH pods and all other pumps work fine. More info is available at [LoopDocs](#).

Follow this [link](#) for more information about checking your **iOS version** or establishing Developer Mode.

Pump

- [Medtronic](#)
- [Omnipod Eros](#)
- [Omnipod Dash](#)
- [Dana\(RS/-i\)](#)
- [Pump Simulator](#)
- *Medtrum (in testing)*

- *Tandem Mobi (in early development)*

CGM

✓ The following CGMs can be used with Trio:

- [Dexcom G5](#)
- [Dexcom G6 / ONE](#)
- [Dexcom G7 / ONE+](#)
- [FreeStyle Libre \(limited models\)](#)
 - US Libre 1 10-day sensors (via transmitters)
 - US Libre 1 14-day sensors (via transmitters)
 - International Libre 1 sensors (via transmitters)
 - European Libre 2 / 2 Plus sensors (via transmitters)
 - European Libre 2 / 2 Plus sensors (direct via Bluetooth)
 - Libre 3 / 3 Plus, all models (requires [LibreLinkUp](#) and continuous internet connection)
- [Eversense E3 / 365](#) (TESTING)
- [Medtronic Enlite](#) (requires a compatible pump)
- [xDrip4iOS](#)
- [\(Nightscout\)](#)
- [\(Glucose Simulator\)](#)

Unsupported FreeStyle Libre sensors

- Non-European FreeStyle Libre 2 (US models not supported)
- FreeStyle Libre Pro
- FreeStyle Libre H

 The table below shows the **compatibility matrix for FreeStyle Libre sensors**.

Model	Region	Duration	Supported?	Via
Libre 1	US	10 days	✓ Yes	Transmitter
Libre 1	US	14 days	✓ Yes	Transmitter
Libre 1	International (non US, non EU)	14 days	✓ Yes	Transmitter

Libre 1	Europe	14 days	 Yes	Transmitter
Libre 2 / 2 Plus	Europe	14 days	 Yes	Bluetooth
Libre 2 (non-European)	US	14 days	 No	N/A
Libre 3	All regions	14 days	 Yes	LibreLinkUp (requires continuous internet connection)
FreeStyle Libre Pro	All regions	14 days	 No	N/A
FreeStyle Libre H	All regions	Varies	 No	N/A

Watches

- [Apple Watch](#)
 - Ultra 1
 - Series 9
 - Series 8
 - Series 7
 - Series 6
 - SE (2nd Gen)
 - Series 5
 - Series 4
- [Garmin Watch](#)