



TrioDocs

Version: 0.4.0

Date: December 12, 2025

Download the latest version on:
<https://triodocs.org>

Share Logs

Sharing Trio Logs

The "Share Logs" feature lets you export the app's logs (from today and yesterday) and opens the iOS share sheet so you can transmit them to support channels, developers, or storage services for troubleshooting and record-keeping.

Check existing solutions

When you run into an issue and want to share your logs, please take a moment to check if it's already been solved by browsing through [Facebook](#) or [Discord](#) first.

Follow the steps below to save and share the logs.

Export Logs

- From within the **Trio** app,
tap  **Settings** > **Support & Communities** > **Share Logs**
- Use the standard **iOS share sheet** to send these files via Mail, cloud-storage apps, or any other destination available in the share sheet.

Once you've exported the logs, you can share them where needed.

Share Logs

Logs are most commonly shared in the following cases:

- **Community troubleshooting:**
Post logs to the Trio [#bug-reports](#) channel on **Discord** or relevant **Facebook** groups when volunteers ask for them to help diagnose an issue.

Overcoming file size limits: Compress Logs

If Discord or Facebook complains about the file size when you upload the logs, **compress them** into an archive before sending.

- **Developer bug reports:** When requested, attach logs to a [GitHub issue](#) so that developers can use them to reproduce, investigate, and resolve the problem.

Compressing Logs

Discord and *Facebook* will **reject** your **file upload** if it is **larger than 25 MB**.

In this case, you first need to compress the log files into an archive before sending.

The following sections describe how to do that on [Windows](#) and [macOS](#).

Compress Logs on Windows

To compress the log files on *Windows*, you can use either a [graphical application](#) or the [terminal](#).

Compress Logs on Windows with a Graphical Application

1. Launch the Windows **File Explorer**

2. Select Log Files

- Navigate to the folder containing the log files.
- Hold `^ Ctrl` and **click** `log.txt`.
- Hold `^ Ctrl` and **click** `log_prev.txt`.

3. Compression

- Right-click one selected file
- Select `Send to` in the contextual menu.
 - Select `Compressed (zipped) folder`.

4. Name Archive

- Type `triologs.zip` and press `Enter`.

Compress Logs on Windows in the Terminal

If you are not comfortable using the **command line in a terminal** on *Windows*, consider using a [graphical application](#) instead; otherwise, keep reading.

Choose your favorite Windows Shell below **PowerShell** or "**Command Prompt (CMD)**", then follow the instructions.

PowerShell

1. Open PowerShell:

- Right-click the `Start` menu
- Select `Windows PowerShell`.

2. Navigate to Log Folder:

```
cd $HOME\Downloads
```

Replace `$HOME/Downloads` with the folder where you downloaded the log files.

3. Run Compression Command:

```
Compress-Archive -Path log.txt, log_prev.txt -DestinationPath triologs.zip
```

where:

- `-Path` is followed by the files to compress.
- `-DestinationPath` is followed by the name of the ZIP archive file.

Compress Logs on macOS

To compress the log files on a *Mac*, you can use either a [graphical application](#) or the [terminal](#).

Compress Logs on macOS with a Graphical Application

1. Launch the **Finder** application
2. **Select Log Files:**
 - Hold `⌘ Cmd` and **click** `log.txt`
 - Hold `⌘ Cmd` and **click** `log_prev.txt`.
3. **Compression :**
 - Right-click the blue highlighted lines.
 - Select **Compress** in the contextual menu.
4. **Name Archive:**
 - An `Archive.zip` file appears.
 - Rename it to `triologs.zip`.

Compress Logs on macOS in the Terminal

If you are not comfortable using the **command line in a terminal** on your Mac, consider using a [graphical application](#) instead; otherwise, keep reading:

1. **Launch Terminal:**
 - Press `⌘ Cmd` + `Space`
 - Type **Terminal**
 - Press `Enter ↵`

2. **Navigate to Downloads Folder:**

```
cd ~/Downloads
```

Replace `~/Downloads` with the folder where you downloaded the log files.

3. **Create Archive:** To create the archive, you can use `zip` or `tar`.

```
zip triologs.zip log.txt log_prev.txt
```

If this method does not work because `zip` is not installed, you can install it like so:

```
xcode-select --install
```

You can also use the *tar gzip alternative* method to compress the log files:

```
tar -czf triologs.tar.gz log.txt log_prev.txt
```