



TrioDocs

Version: 0.4.0

Date: June 15, 2026

Download the latest version on:  
<https://triodocs.org>

# Nightscout

# Nightscout

*Nightscout* is an excellent tool for remotely viewing *Trio* actions. *Nightscout* can act as a stand-alone tool or be integrated with *Trio*. (*Nightscout* integrates with many open-source hybrid closed loop systems such as *Loop*, OpenAPS, AndroidAPS and *Trio*.)

When integrated with *Trio*, *Nightscout* provides monitoring of *Trio* activities such as viewing history of glucose, carbs, boluses, temp basals and overrides; troubleshooting *Trio* errors; and provides extensive reports for analyzing data trends and patterns. These reports assist when *Trio* Therapy Settings need to be adjusted.

For caregivers, *Nightscout* enables remote monitoring and even the ability to issue remote commands through *Nightscout* when both *Trio* and *Nightscout* are properly configured. If you add *LoopFollow*, those same remote commands can be issued from the caregiver's *LoopFollow* phone with better connectivity and feedback than using *Nightscout* alone.

For more information about *Nightscout*, please refer to [Nightscout Documentation](#).

For more information about using remote control with a combination of *Nightscout* and *Trio*, please refer to the [LoopFollowDocs](#).