



TrioDocs

Version: 0.4.0

Date: August 13, 2026

Download the latest version on:
<https://triodocs.org>

Mac Build Overview

Build Trio with Mac - Overview

Congratulations on choosing to build Trio using your Mac computer. This documentation provides all the necessary instructions to download a local copy of the Trio repository using the loopandlearn `TrioBuildSelectScript`. It will then guide you on using a compatible computer with the Xcode application to build and deploy the Trio application to a compatible device.

Requirements:

- [Compatible Phone ↗](#)
- [Compatible CGM ↗](#)
- [Compatible Insulin Pump ↗](#)
- [Apple Developer Account ↗](#)
- [Compatible Computer ↗](#)
- [Xcode ↗](#)

Getting Ready to Build

The first step in getting ready to build Trio using your Mac computer is to review the above requirements. These are considered the minimum requirements for Building Trio using your Mac. Don't stress if you are unable to meet the above requirements. Other options are available for building Trio, and I recommend you start with [Build Trio with Browser ↗](#).

The second step is to read over and become familiar with the build instructions in their entirety. This includes familiarizing yourself with common [Build Errors ↗](#). Most mistakes you can make while building Trio have already been made, and our dedicated community has documented necessary fixes.

When you are ready to proceed, take your time, read carefully, and complete the tasks on each page. Read the top three boxes on each page; icons for those boxes are displayed below for reference.

Overview

Time Estimate

- 5 minutes, if your computer meets minimum requirements.
- 30-60 minutes, if you need to install macOS updates.

Summary

To build the Trio app on a Mac, your computer, iPhone, and Xcode must have compatible versions. If you are buying a Mac to use the build with the Mac method, choose one that can be updated to the Tahoe (macOS 26) operating system and has at least 256 GB (512 GB is better). The Build with Browser method works on any computer or tablet.

? FAQs

- **Do I need a Mac?** No, you can [Build Trio with Browser](#) on any computer.
- **How often do I need to use the computer to build with the Mac method?** Computer access is required when:
 - Initially installing the Trio app.
 - Trio will expire annually for a paid account or weekly for a free account.
 - Updating to a newer Trio release.
 - You do not need access to an Apple computer to update your phone's iOS, troubleshoot, or change Trio settings
- **What if my computer doesn't meet the minimum requirements** If your 'macOS' does not meet the minimum standards as your iOS version dictates, you will need to complete a software update on your computer to continue. Suppose you cannot update your computer to the required 'macOS' due to hardware limitations. In that case, your options are as follows:
 - Build Trio with GitHub using the Browser Build method.
 - Borrow an up-to-date Mac computer.
 - Purchase an up-to-date Mac computer.

? What if I get stuck?

- First, take some time to read over the instructions to ensure you have understood all the steps.
- Compare your screen's display with the images in the documents. Is there something different? Does yours have an error message? If so, does it guide you to the problem and solution?
- If you are still stuck, you can seek assistance via [Facebook](#) or [Discord](#).